

BRIGHTON VOLUNTEER AMBULANCE

VENTILATOR/BAG CHECKLIST

<p>POWER CABLE (1) _____</p> <p>HIGH PRESSURE OXYGEN CABLE (1) _____</p> <p>SUCTION TUBING (1) _____</p> <p>14 FRENCH SUCTION CATHETER (1) _____</p> <p>IN-LINE ETCO2 (2) _____</p>	<p>FRONT POUCH</p>	<p>PROVIDER NAME</p> <p>_____</p> <p>DATE</p> <p>_____</p> <p>***BOLD/ITALICIZED ITEMS REQUIRE EXPIRATION DATES***</p>
<p>OXYGEN REGULATOR (DISS PORT) (1) _____</p> <p>OXYGEN D TANK (MIN 2000 PSI) (1) _____</p> <p>ADULT VENTILATOR CIRCUIT (2) _____</p>	<p>CENTER POUCH</p>	
<p>PEEP VALVE (1) _____</p> <p>PEDIATRIC BVM (1) _____</p>	<p>RIGHT SIDED POUCH</p>	
<p>PEDIATRIC VENTILATOR CIRCUIT (1) _____</p>	<p>REAR POUCH</p>	
<p>HIGH PRESSURE OXYGEN SWIZEL PORT (1) _____</p>	<p>LEFT SIDED POUCH</p>	
<p>DEVICE CHECK (TURN ON/OFF) _____</p> <p>CARRYING SHOULDER STRAP _____</p> <p>BATTERY FULLY CHARGED _____</p> <p>PLUGGED IN/CHARGING _____</p>	<p>AHP 300 VENT</p>	